

Fibulin-1

Catalog No: #AP83110



Package Size: #AP83110-1 50ug #AP83110-2 100ug #AP83110-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fibulin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FIBL-1,Basement-membrane protein 90,BM-90
Accession No.	Uniprot ID:Q08879
Uniprot	Q08879
GeneID	14114;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fbln1

Note: This product is for in vitro research use only