

Ferritin light chain 2

Catalog No: #AP83112



Package Size: #AP83112-1 50ug #AP83112-2 100ug #AP83112-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ferritin light chain 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Ferritin L subunit 2,Ferritin subunit LG,Ftl-2
Accession No.	Uniprot ID:P49945
Uniprot	P49945
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ftl2

Note: This product is for in vitro research use only