

## Alpha-ketoglutarate-dependent dioxygenase FTO

Catalog No: #AP83114



Package Size: #AP83114-1 50ug #AP83114-2 100ug #AP83114-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Alpha-ketoglutarate-dependent dioxygenase FTO
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Fat mass and obesity-associated protein
Accession No.	Uniprot ID:Q2A121
Uniprot	Q2A121
GeneID	291905;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Fto

Note: This product is for in vitro research use only