

Keratinocyte growth factor

Catalog No: #AP83115



Package Size: #AP83115-1 50ug #AP83115-2 100ug #AP83115-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Keratinocyte growth factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	KGF,Fibroblast growth factor 7,FGF-7,Heparin-binding growth factor 7,HBGF-7,KGF
Accession No.	Uniprot ID:P79150
Uniprot	P79150
GeneID	403915;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FGF7

Note: This product is for in vitro research use only