

Four and a half LIM domains protein 2

Catalog No: #AP83116



Package Size: #AP83116-1 50ug #AP83116-2 100ug #AP83116-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Four and a half LIM domains protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FHL-2,Skeletal muscle LIM-protein 3,SLIM-3
Accession No.	Uniprot ID:O70433
Uniprot	O70433
GeneID	14200;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fhl2

Note: This product is for in vitro research use only