

DAPK1 (Phospho-Ser289) Antibody

Catalog No: #12496



Package Size: #12496-1 50ul #12496-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DAPK1 (Phospho-Ser289) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms
Specificity	DAPK1 (Phospho-Ser289) Antibody detects endogenous levels of DAPK1 only when phosphorylated at Ser289
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human DAPK1 (Phospho-Ser289)
Target Name	DAPK1
Modification	Phospho
Other Names	DAPK1, DAP kinase 1, DAPK
Accession No.	Swiss-Prot#: P53355NCBI Gene ID: 1612
Uniprot	P53355
GeneID	1612;
Target Species	human
Calculated MW	160kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images



Western Blot analysis of lysates of 1. HeLa cell and 2. HeLa cell treated with 100ng/mL LPS for 30min, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only