

Glutathione S-transferase A4

Catalog No: #AP83121



Package Size: #AP83121-1 50ug #AP83121-2 100ug #AP83121-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Glutathione S-transferase A4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	GST class-alpha member 4, Glutathione S-transferase alpha-4
Accession No.	Uniprot ID:Q5E9G0
Uniprot	Q5E9G0
GeneID	533917;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GSTA4

Note: This product is for in vitro research use only