

Endothelin-2

Catalog No: #AP83122



Package Size: #AP83122-1 50ug #AP83122-2 100ug #AP83122-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Endothelin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	ET-2,Preproendothelin-2,PPET2
Accession No.	Uniprot ID:P12064
Uniprot	P12064
GeneID	403508;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EDN2

Note: This product is for in vitro research use only