

# Fibroblast growth factor 7

Catalog No: #AP83128



Package Size: #AP83128-1 50ug #AP83128-2 100ug #AP83128-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Fibroblast growth factor 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FGF-7,Heparin-binding growth factor 7,HBGF-7,Keratinocyte growth factor,KGF
Accession No.	Uniprot ID:P21781
Uniprot	P21781
GeneID	2252;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:FGF7

Note: This product is for in vitro research use only