

# Glycoprotein hormones alpha chain

Catalog No: #AP83130

Package Size: #AP83130-1 50ug #AP83130-2 100ug #AP83130-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Glycoprotein hormones alpha chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Anterior pituitary glycoprotein hormones common subunit alpha,Follicle-stimulating hormone alpha chain,FSH-alpha,Follitropin alpha chain,Luteinizing hormone alpha chain,LSH-alpha,Lutropin alpha chain,Thyroid-stimulating hormone alpha chain,TSH-alpha,Thyrotropin alpha chain
Accession No.	Uniprot ID:Q9XSW8
Uniprot	Q9XSW8
GeneID	403483;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CGA

Note: This product is for in vitro research use only