

High affinity immunoglobulin epsilon receptor subunit alpha



Catalog No: #AP83132

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP83132-1 50ug #AP83132-2 100ug #AP83132-3 1mg

Description

Product Name	High affinity immunoglobulin epsilon receptor subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Fc-epsilon RI-alpha,FcERI,IgE Fc receptor subunit alpha,FCE1A
Accession No.	Uniprot ID:P12319
Uniprot	P12319
GeneID	2205;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FCER1A

Note: This product is for in vitro research use only