

Fibroblast growth factor 8

Catalog No: #AP83134



Package Size: #AP83134-1 50ug #AP83134-2 100ug #AP83134-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fibroblast growth factor 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	FGF-8,Heparin-binding growth factor 8,HBGF-8
Accession No.	Uniprot ID:Q90722
Uniprot	Q90722
GeneID	396313;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FGF8

Note: This product is for in vitro research use only