

## Four and a half LIM domains protein 1

Catalog No: #AP83139



Package Size: #AP83139-1 50ug #AP83139-2 100ug #AP83139-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Four and a half LIM domains protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FHL-1,KyoT,RBP-associated molecule 14-1, RAM14-1,Skeletal muscle LIM-protein 1,SLIM
Accession No.	Uniprot ID:P97447
Uniprot	P97447
GeneID	14199;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Fhl1

Note: This product is for in vitro research use only