

## Endoplasmin

Catalog No: #AP83142



Package Size: #AP83142-1 50ug #AP83142-2 100ug #AP83142-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Endoplasmin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	94 kDa glucose-regulated protein,GRP-94,Heat shock protein 90 kDa beta member 1,GRP94,TRA1
Accession No.	Uniprot ID:Q95M18
Uniprot	Q95M18
GeneID	282646;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HSP90B1

Note: This product is for in vitro research use only