

Squalene monooxygenase

Catalog No: #AP83144



Package Size: #AP83144-1 50ug #AP83144-2 100ug #AP83144-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Squalene monooxygenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Squalene epoxidase,SE,Erg1
Accession No.	Uniprot ID:P52020
Uniprot	P52020
GeneID	29230;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sqle

Note: This product is for in vitro research use only