

Coagulation factor XI

Catalog No: #AP83147



Package Size: #AP83147-1 50ug #AP83147-2 100ug #AP83147-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Coagulation factor XI
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	FXI, Plasma thromboplastin antecedent, PTA
Accession No.	Uniprot ID:Q5NTB3
Uniprot	Q5NTB3
GeneID	407998;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: F11

Note: This product is for in vitro research use only