

Follicle-stimulating hormone receptor

Catalog No: #AP83148



Package Size: #AP83148-1 50ug #AP83148-2 100ug #AP83148-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Follicle-stimulating hormone receptor |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Pig |
| Other Names | FSH-R,Follitropin receptor |
| Accession No. | Uniprot ID:P49059 |
| Uniprot | P49059 |
| GeneID | 397679; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FSHR

Note: This product is for in vitro research use only