

Alpha-(1,3)-fucosyltransferase 6

Catalog No: #AP83150



Package Size: #AP83150-1 50ug #AP83150-2 100ug #AP83150-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-(1,3)-fucosyltransferase 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Fucosyltransferase 6,Fucosyltransferase VI,Fuc-TVI,Galactoside 3-L-fucosyltransferase,FCT3A
Accession No.	Uniprot ID:P51993
Uniprot	P51993
GeneID	2528;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FUT6

Note: This product is for in vitro research use only