

Glutaryl-CoA dehydrogenase, mitochondrial

Catalog No: #AP83151



Package Size: #AP83151-1 50ug #AP83151-2 100ug #AP83151-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glutaryl-CoA dehydrogenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	GCD
Accession No.	Uniprot ID:Q2KHZ9
Uniprot	Q2KHZ9
GeneID	506310;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GCDH

Note: This product is for in vitro research use only