

Glucocorticoid receptor

Catalog No: #AP83152



Package Size: #AP83152-1 50ug #AP83152-2 100ug #AP83152-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Glucocorticoid receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	GR,Nuclear receptor subfamily 3 group C member 1,Grl
Accession No.	Uniprot ID:P06536
Uniprot	P06536
GeneID	24413;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nr3c1

Note: This product is for in vitro research use only