

# FAS-associated factor 2

Catalog No: #AP83154



Package Size: #AP83154-1 50ug #AP83154-2 100ug #AP83154-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FAS-associated factor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	UBX domain-containing protein 8,Kiaa0887,Ubx8
Accession No.	Uniprot ID:Q3TDN2
Uniprot	Q3TDN2
GeneID	76577;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Faf2

Note: This product is for in vitro research use only