

# Fibroblast growth factor 21

Catalog No: #AP83155



Package Size: #AP83155-1 50ug #AP83155-2 100ug #AP83155-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Fibroblast growth factor 21
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FGF-21,UNQ3115/PRO10196
Accession No.	Uniprot ID:Q9NSA1
Uniprot	Q9NSA1
GeneID	26291;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:FGF21

Note: This product is for in vitro research use only