

Forkhead box protein C2

Catalog No: #AP83158



Package Size: #AP83158-1 50ug #AP83158-2 100ug #AP83158-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Forkhead box protein C2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Brain factor 3,BF-3,HFH-BF-3,Hfhbf3
Accession No.	Uniprot ID:Q63246
Uniprot	Q63246
GeneID	171356;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Foxc2

Note: This product is for in vitro research use only