

Alpha-2-HS-glycoprotein

Catalog No: #AP83160



Package Size: #AP83160-1 50ug #AP83160-2 100ug #AP83160-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alpha-2-HS-glycoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Countertrypin,Fetuin-A,Fetua
Accession No.	Uniprot ID:P29699
Uniprot	P29699
GeneID	11625;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ahsg

Note: This product is for in vitro research use only