

## Endothelial cell-specific molecule 1

Catalog No: #AP83163

Package Size: #AP83163-1 50ug #AP83163-2 100ug #AP83163-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Endothelial cell-specific molecule 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	ESM-1,PG25
Accession No.	Uniprot ID:P97682
Uniprot	P97682
GeneID	64536;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Esm1

Note: This product is for in vitro research use only