

## Furin

Catalog No: #AP83166



Package Size: #AP83166-1 50ug #AP83166-2 100ug #AP83166-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Furin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Dibasic-processing enzyme, Paired basic amino acid residue-cleaving enzyme, PACE, Trans Golgi network protease furin, FUR, PACE
Accession No.	Uniprot ID: Q28193
Uniprot	Q28193
GeneID	281374;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: FURIN

Note: This product is for in vitro research use only