

Epithelial membrane protein 1

Catalog No: #AP83169



Package Size: #AP83169-1 50ug #AP83169-2 100ug #AP83169-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Epithelial membrane protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	EMP-1, Tumor-associated membrane protein
Accession No.	Uniprot ID:P54848
Uniprot	P54848
GeneID	25314;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Emp1

Note: This product is for in vitro research use only