

## Ephrin type-A receptor 3

Catalog No: #AP83170



Package Size: #AP83170-1 50ug #AP83170-2 100ug #AP83170-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ephrin type-A receptor 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	EPH-like kinase 4,EK4,HEK,Human embryo kinase,Tyrosine-protein kinase TYRO4,Tyrosine-protein kinase receptor ETK1,Eph-like tyrosine kinase 1,ETK,ETK1,HEK,TYRO4
Accession No.	Uniprot ID:P29320
Uniprot	P29320
GeneID	2042;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:EPHA3

Note: This product is for in vitro research use only