

Ewing's tumor-associated antigen 1 homolog

Catalog No: #AP83171



Package Size: #AP83171-1 50ug #AP83171-2 100ug #AP83171-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ewing's tumor-associated antigen 1 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:Q08DI1
Uniprot	Q08DI1
GeneID	282224;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ETAA1

Note: This product is for in vitro research use only