

Free fatty acid receptor 2

Catalog No: #AP83172



Package Size: #AP83172-1 50ug #AP83172-2 100ug #AP83172-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Free fatty acid receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	G-protein coupled receptor 43,Gpr43
Accession No.	Uniprot ID:Q76E16
Uniprot	Q76E16
GeneID	292794;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ffar2

Note: This product is for in vitro research use only