

Electron transfer flavoprotein subunit beta

Catalog No: #AP83174



Package Size: #AP83174-1 50ug #AP83174-2 100ug #AP83174-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Electron transfer flavoprotein subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Beta-ETF
Accession No.	Uniprot ID:Q6UAQ8
Uniprot	Q6UAQ8
GeneID	396614;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ETFB

Note: This product is for in vitro research use only