

Coagulation factor XI

Catalog No: #AP83176



Package Size: #AP83176-1 50ug #AP83176-2 100ug #AP83176-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Coagulation factor XI
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FXI, Plasma thromboplastin antecedent, PTA
Accession No.	Uniprot ID:Q91Y47
Uniprot	Q91Y47
GeneID	109821;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: F11

Note: This product is for in vitro research use only