

# Glutamate receptor 1

Catalog No: #AP83181



Package Size: #AP83181-1 50ug #AP83181-2 100ug #AP83181-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Glutamate receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	GluR-1,AMPA-selective glutamate receptor 1,GluR-A,GluR-K1,Glutamate receptor ionotropic, AMPA 1,GluA1,Glur1
Accession No.	Uniprot ID:P19490
Uniprot	P19490
GeneID	50592;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Gria1

Note: This product is for in vitro research use only