

IgGFc-binding protein

Catalog No: #AP83183



Package Size: #AP83183-1 50ug #AP83183-2 100ug #AP83183-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	IgGFc-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Fcgamma-binding protein antigen,FcgammaBP
Accession No.	Uniprot ID:Q9Y6R7
Uniprot	Q9Y6R7
GeneID	8857;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FCGBP

Note: This product is for in vitro research use only