

Fibroblast growth factor 4

Catalog No: #AP83185



Package Size: #AP83185-1 50ug #AP83185-2 100ug #AP83185-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibroblast growth factor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FGF-4,Heparin secretory-transforming protein 1,HST,Heparin-binding growth factor 4,HBGF-4,Transforming protein KS3,HST,HSTF1,KS3
Accession No.	Uniprot ID:P08620
Uniprot	P08620
GeneID	2249;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FGF4

Note: This product is for in vitro research use only