

# EGF-like repeat and discoidin I-like domain-containing protein 3



Catalog No: #AP83186

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP83186-1 50ug #AP83186-2 100ug #AP83186-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | EGF-like repeat and discoidin I-like domain-containing protein 3                              |
| Brief Description  | Recombinant Protein   |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE  |
| Species Reactivity | Hu  |
| Other Names        | Developmentally-regulated endothelial cell locus 1 protein,Integrin-binding protein DEL1,DEL1 |
| Accession No.      | Uniprot ID:O43854   |
| Uniprot            | O43854  |
| GeneID             | 10085;  |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)  |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:EDIL3

Note: This product is for in vitro research use only