

Gamma-interferon-inducible lysosomal thiol reductase

Catalog No: #AP83187



Package Size: #AP83187-1 50ug #AP83187-2 100ug #AP83187-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Gamma-interferon-inducible lysosomal thiol reductase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Interferon gamma inducible protein 30
Accession No.	Uniprot ID:Q499T2
Uniprot	Q499T2
GeneID	290644;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ifi30

Note: This product is for in vitro research use only