

## F-box only protein 32

Catalog No: #AP83190



Package Size: #AP83190-1 50ug #AP83190-2 100ug #AP83190-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	F-box only protein 32
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Atrogin-1, Muscle atrophy F-box protein, MAFbx
Accession No.	Uniprot ID:Q9CPU7
Uniprot	Q9CPU7
GeneID	67731;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Fbxo32

Note: This product is for in vitro research use only