

Alpha-fetoprotein

Catalog No: #AP83193



Package Size: #AP83193-1 50ug #AP83193-2 100ug #AP83193-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alpha-fetoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Alpha-1-fetoprotein,Alpha-fetoglobulin
Accession No.	Uniprot ID:Q8MJU5
Uniprot	Q8MJU5
GeneID	403551;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AFP

Note: This product is for in vitro research use only