

Fibroblast growth factor-binding protein 2

Catalog No: #AP83194



Package Size: #AP83194-1 50ug #AP83194-2 100ug #AP83194-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Fibroblast growth factor-binding protein 2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | FGF-BP2,37 kDa killer-specific secretory protein,Ksp37,HBp17-related protein,HBp17-RP,KSP37,UNQ425/PRO1065 |
| Accession No. | Uniprot ID:Q9BYJ0 |
| Uniprot | Q9BYJ0 |
| GeneID | 83888; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FGFBP2

Note: This product is for in vitro research use only