

Double-stranded RNA-specific adenosine deaminase

Catalog No: #AP83197



Package Size: #AP83197-1 50ug #AP83197-2 100ug #AP83197-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Double-stranded RNA-specific adenosine deaminase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	DRADA,Dsrad
Accession No.	Uniprot ID:P55266
Uniprot	P55266
GeneID	81635;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Adar

Note: This product is for in vitro research use only