

Dynamin-1

Catalog No: #AP83198



Package Size: #AP83198-1 50ug #AP83198-2 100ug #AP83198-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dynamin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	B-dynamin,D100,Dynamin, brain,Dnm
Accession No.	Uniprot ID:P21575
Uniprot	P21575
GeneID	140694;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dnm1

Note: This product is for in vitro research use only