

Glycosylation-dependent cell adhesion molecule 1

Catalog No: #AP83199



Package Size: #AP83199-1 50ug #AP83199-2 100ug #AP83199-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glycosylation-dependent cell adhesion molecule 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	GlyCAM-1, Endothelial ligand FOR L-selectin, MC26, SGP50, Sulfated 50 kDa glycoprotein
Accession No.	Uniprot ID:Q02596
Uniprot	Q02596
GeneID	14663;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Glycam1

Note: This product is for in vitro research use only