

FGFR3 (Phospho-Tyr760) Antibody

Catalog No: #12504



Package Size: #12504-1 50ul #12504-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FGFR3 (Phospho-Tyr760) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	FGFR3 (Phospho-Tyr760) Antibody detects endogenous levels of FGFR3 only when phosphorylated at Tyr760
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human FGFR3 (Phospho-Tyr760)
Target Name	FGFR3
Modification	Phospho
Other Names	FGFR3, CD333, CD333 antigen, ACH, FGF receptor 3, HSGFR3EX, Hydroxyaryl-protein kinase, JTK4, FGFR-3, CEK2, Tyrosine kinase JTK4
Accession No.	Swiss-Prot#: P22607 NCBI Gene ID: 2261
Uniprot	P22607
GeneID	2261;
Target Species	human
Calculated MW	95kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images



Western Blot analysis of lysates of SY-SY5Y cell, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only