

Follistatin-related protein 3

Catalog No: #AP83212



Package Size: #AP83212-1 50ug #AP83212-2 100ug #AP83212-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Follistatin-related protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Follistatin-like protein 3,Follistatin-related gene protein,Flrg
Accession No.	Uniprot ID:Q9EQC7
Uniprot	Q9EQC7
GeneID	83554;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fstl3

Note: This product is for in vitro research use only