

## Glyceraldehyde-3-phosphate dehydrogenase

Catalog No: #AP83213



Package Size: #AP83213-1 50ug #AP83213-2 100ug #AP83213-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Glyceraldehyde-3-phosphate dehydrogenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	GAPDH,Peptidyl-cysteine S-nitrosylase GAPDH,Gapd
Accession No.	Uniprot ID:P16858
Uniprot	P16858
GeneID	100042025;14433;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Gapdh

Note: This product is for in vitro research use only