

Proheparin-binding EGF-like growth factor

Catalog No: #AP83216



Package Size: #AP83216-1 50ug #AP83216-2 100ug #AP83216-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Proheparin-binding EGF-like growth factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Dtr,Hegfl
Accession No.	Uniprot ID:Q06186
Uniprot	Q06186
GeneID	15200;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hbegf

Note: This product is for in vitro research use only