

Alpha-enolase

Catalog No: #AP83217



Package Size: #AP83217-1 50ug #AP83217-2 100ug #AP83217-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-enolase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	2-phospho-D-glycerate hydro-lyase,Enolase 1,Non-neural enolase,NNE,Eno-1
Accession No.	Uniprot ID:P04764
Uniprot	P04764
GeneID	24333;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Eno1

Note: This product is for in vitro research use only