

Ecto-NOX disulfide-thiol exchanger 2

Catalog No: #AP83218



Package Size: #AP83218-1 50ug #AP83218-2 100ug #AP83218-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ecto-NOX disulfide-thiol exchanger 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	APK1 antigen,Cytosolic ovarian carcinoma antigen 1,Tumor-associated hydroquinone oxidase,tNOX,Cova1
Accession No.	Uniprot ID:Q8R0Z2
Uniprot	Q8R0Z2
GeneID	209224;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Enox2

Note: This product is for in vitro research use only