

Proto-oncogene c-Fos

Catalog No: #AP83220



Package Size: #AP83220-1 50ug #AP83220-2 100ug #AP83220-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proto-oncogene c-Fos
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cellular oncogene fos,G0/G1 switch regulatory protein 7,G0S7
Accession No.	Uniprot ID:P01100
Uniprot	P01100
GeneID	2353;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FOS

Note: This product is for in vitro research use only